

Immunology Consultants Laboratory, Inc.

15862 SW 72nd Avenue, #150
Portland, OR 97224
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

- PRODUCT:** Agarose Immobilized Rabbit anti-C-MYC
- CATALOG NUMBER:** RMYC-145D-Z
- LOT NUMBER:**
- QUANTITY:** 100 µg per 200 µl gel.
- PHYSICAL STATE:** Liquid Suspension – Total volume 400 µl
- BUFFER:** Phosphate buffered saline, pH 7.2
- PRESERVATIVE:** 0.1% Sodium Azide
- STORAGE:** 2°- 4° C
- EXPIRATION:** 12 months after shipment date.
- PROCEDURE:** Rabbits were immunized with C-MYC peptide conjugated with KLH. After multiple immunizations in Freund's adjuvant, serum was collected. Antibodies were immunoaffinity purified using the peptide immobilized on a solid phase then coupled to agarose beads using aldehyde chemistry.
- APPLICATIONS:** Immunoprecipitation: 15 – 25 ul of gel slurry per 0.2 – 1.0 mg of lysate or extract.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

CERTIFICATE OF ANALYSIS VERIFICATION:
PRODUCT OF THE UNITED STATES OF AMERICA

DATE: _____

RMYC145DZ