

Certificate of Analysis

Product Code:	MA-FD1
Lot Number:	32049
Antibody:	FD1 anti Fel d 1 mAb
Immunogen:	Fel d 1
Isotype	Mouse IgG
Purification	Produced in ascites and purified by chromatography using Protein G. Single heavy and light chain bands on SDS-PAGE.
Concentration:	4.57 mg/ml , based on A280 for IgG (1.42= 1mg/ml) .
Quantity:	200µl
Application:	Use at 1:1000 in ELISA. Final working dilutions for other applications must be determined by end user.
Presentation:	In phosphate buffered saline, pH 7.2. 0.22µm filtered, preservative free.
Storage/Handling:	Maintain at 4°C for up to 12 months.

Indoor Biotechnologies, Inc.
1216 Harris Street,
Charlottesville
Virginia, 22903
United States

Tel: (434) 984-2304
Fax: (434) 984-2709
mail@inbio.com

Indoor Biotechnologies Ltd
The Old Brewery
38, High Street , Warminster
Wiltshire, BA12 9AF
United Kingdom

Tel: 44 (0)5601 153 291
Fax: 44 (0)1985 218 300
info@indoorbiotech.co.uk