



## Recombinant Cyp c 1

**Product Code:** RE-CC1-1

**Gel Image:**

**Allergen:** rCyp c 1 (*Cyprinus carpio* Allergen 1.0101)

**Lot No:** XXXX

**Source:** *E. coli*

**Mol. Wt:** 12 kD

**Purification:** Multi-step chromatographic purification.  
Purity > 95 % by SDS-PAGE.

**Concentration:** See Product Insert

**Formulation:** Preservative-free and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** See Product Insert

**Notes:** (1) rCyp c 1 has a C-terminal 6xHis-tag.



**For research and commercial use in vitro: not for human in vivo or therapeutic use.**

**An InBio® product. Made in the USA**

Indoor Biotechnologies, Inc.  
700 Harris Street  
Charlottesville, VA 22903, USA  
www.inbio.com

InBio® US: mail@inbio.com  
InBio® UK: info@indoorbiotech.co.uk