



DATA SHEET

HUMAN PMN ELASTASE

Catalog no.:	A 6841AG.1 / A 6841AG.2
Protein:	Human PMN Elastase, purified
Reference no.:	LocusID 4879
Source:	Human Sputum
Purity:	> 95% (SDS-PAGE)
Contents:	100 µg/1 mg (lyophilized) Resuspend in 0,05 M NaOAc 0,1 M NaCl pH 5,00 for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD.DAT\DATEN\BLAETTER\PROTEINE\A 1005AG-proBNP1-76rec.DOC