



Product:	Mouse Factor X
Catalog #:	IMFX
Lot Number:	Sample
Concentration:	6.1mg/ml
Specific Activity:	110 U/mg
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 11.6$
Molecular Weight:	58900
Buffer:	50% Glycerol/ H <sub>2</sub> O (vol/vol)
Recommended Storage:	-20 ° C

**PRODUCT/CAUTION NOTES:**

Products contained within this shipment are of human blood/plasma origin. Although the starting material was tested prior to initiation of the manufacturing process, and was found negative or nonreactive for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II, extreme caution should be used when handling this material as there is a margin of error in all tests.

This product is to be used for IN VITRO research purposes only, and is not to be used for clinical, diagnostic or human use!

\*Determined by Factor X clotting assay. One unit is equivalent to the Factor X activity in one milliliter of normal plasma.