

## Product Specifications

Product:	UltraPolymer Donkey anti-Sheep IgG (H&L) – HRP
Description:	UltraPolymer Donkey anti-Sheep IgG (H&L) conjugated to HRP, affinity purified, min x w/human, mouse or rabbit IgG
Part Number:	DASHRP-050
Concentration:	n/a
Amount:	50 ml
Conjugate:	Horseradish Peroxidase
Form:	Liquid
Purification:	Affinity purified using solid phase Sheep IgG
Purity:	Affinity purified antibody is > 95% based on SDS-PAGE
Host:	Donkey
Immunogen:	Purified Sheep IgG, whole molecule
Buffer:	PBS, 1% BSA & 0.1% proclin 150
Storage:	2-8 °C
Shelf Life:	1 year from date of receipt. Prepare working dilution prior to use and then discard.
Specificity:	Based on IEP, this immunopolymer antibody reacts with heavy ( $\gamma$ ) chains on sheep IgG, light chains on all sheep immunoglobulins
Cross Reactivity:	Based on IEP, no cross-reactivity is observed to non-immunoglobulin sheep serum proteins, IgG from human, mouse or rabbit
Country of Origin:	Donkey serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.
Applications:	For use as a standard or control in most immunoassay formats.
Disclaimer:	For <i>in vitro</i> Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. <b>Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.</b>