

## Product Specifications

Product:	Donkey anti-Mouse IgG (H&L), F(ab)'2 Fragment
Description:	Donkey anti-Mouse IgG (H&L) - F(ab)'2 Fragment, DyLight®633 Conjugate, min x w/bovine, chicken, goat, guinea pig, hamster, horse, human, rabbit, rat or sheep IgG
Part Number:	DkxMu-003-L633NHSX
Concentration:	1 mg/ml (E 1% at 280 nm = 13.0)
Amount:	0.5 mg
Conjugate:	DyLight® 633 (Ex = 638 nm; Em = 658 nm)
Form:	Lyophilized
Purification:	Affinity purified using solid phase Mouse IgG
Purity:	Affinity purified antibody is > 95% based on SDS-PAGE
Host:	Donkey
Immunogen:	Purified Mouse IgG, whole molecule
Buffer:	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free
Preservative:	0.05% (w/v) Sodium Azide
Reconstitution:	Rehydrate with 0.55 ml of deionized water and let stand 30 minutes at room temperature to dissolve. Centrifuge to remove any particulates. Prepare fresh working dilution daily.
Storage:	Store freeze-dried powder at 2-8 °C.
Shelf Life:	Product is stable for up to 4 weeks at 2-8°C after rehydration. For extended storage after rehydration, add an equal volume of glycerol and store at -20°C.
Specificity:	Based on IEP, this antibody reacts with: <ul style="list-style-type: none"><li>• heavy (γ) chains on mouse IgG</li><li>• light chains on all mouse immunoglobulins</li></ul>
Cross Reactivity:	Based on IEP, no reactivity is observed to: <ul style="list-style-type: none"><li>• non-immunoglobulin mouse serum immunoglobulins</li><li>• IgG from bovine, chicken, goat, guinea pig, hamster, horse, human, rabbit, rat or sheep</li></ul>
Country of Origin:	Donkey serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.
Applications:	This conjugate is suitable for immunomicroscopy and flow cytometry. The optimal working dilution should be determined by the investigator. Suggested starting

dilution(s):

- 1:20 - 1:2,000 for most applications

Disclaimer:

For *in vitro* Laboratory Use Only. Not intended for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.

**Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.**

Trademark:

DyLight® is a trademark of Thermo Fisher Scientific, Inc. and its subsidiaries.