

## Product Specifications

Product:	UltraPolymer Goat anti-Mouse IgG (H&L) – HRP
Description:	UltraPolymer Goat anti-Mouse IgG (H&L) conjugated to HRP, affinity purified, min x w/bovine, horse, human, pig or rabbit serum proteins
Part Number:	GAMHRP-050
Concentration:	n/a
Amount:	50 ml
Conjugate:	Horseradish Peroxidase
Form:	Liquid
Purification:	Affinity purified using solid phase Mouse IgG
Purity:	Affinity purified antibody is > 95% based on SDS-PAGE
Host:	Goat
Immunogen:	Purified Mouse IgG, whole molecule
Buffer:	PBS & 0.1% proclin 150
Storage:	2-8 °C
Specificity:	Based on IEP, this immunopolymer antibody reacts with heavy ( $\gamma$ ) chains on mouse IgG, light chains on all mouse immunoglobulins
Cross Reactivity:	Based on IEP, no cross-reactivity is observed to non-immunoglobulin mouse serum proteins, serum proteins from bovine, horse, human, pig or rabbit
Country of Origin:	Goat serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.
Applications:	For use as a standard or control in most immunoassay formats.
Disclaimer:	For <i>in vitro</i> Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. <b>Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.</b>