

Product Specifications

Product:	Goat anti-Guinea Pig IgG (H&L) - Affinity Pure	
Description:	Goat anti-Guinea Pig IgG (H&L) - Affinity Pure, DyLight®594 Conjugate	
Part Number:	GtxGp-003-D594NHSX	
Concentration:	1.0 mg/ml (E 1% at 280 nm = 13.0)	
Amount:	1.0 mg	
Conjugate:	DyLight® 594 (Ex = 593 nm; Em = 618 nm)	
Form:	Lyophilized	
Purification:	Affinity purified using solid phase Guinea Pig IgG	
Purity:	Affinity purified antibody is > 95% based on SDS-PAGE	
Host:	Goat	
Immunogen:	Purified Guinea Pig IgG, whole molecule	
Buffer:	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free	
Preservative:	0.05% (w/v) Sodium Azide	
Reconstitution:	Rehydrate with 1.1 ml of deionized water and let stand 30 minutes at room temperature to dissolve. Centrifuge to remove any particulates. Prepare fresh working dilution daily.	
Storage:	Store freeze-dried powder at 2-8 °C.	
Shelf Life:	Product is stable for up to 4 weeks at 2-8°C after rehydration. For extended storage after rehydration, add an equal volume of glycerol and store at -20°C.	
Specificity:	Based on IEP, this antibody reacts with: <ul style="list-style-type: none">● heavy (γ) chains on guinea pig IgG● light chains on all guinea pig immunoglobulins	
Cross Reactivity:	Based on IEP, no reactivity is observed to: <ul style="list-style-type: none">● non-immunoglobulin guinea pig serum immunoglobulins	
Country of Origin:	Goat serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.	
Applications:	This conjugate is suitable for immunomicroscopy and flow cytometry. The optimal working dilution should be determined by the investigator. Suggested starting dilution(s): <ul style="list-style-type: none">● 1:20 - 1:2,000 for most applications	
Disclaimer:	For <i>in vitro</i> Laboratory Use Only. Not intended for diagnostic or therapeutic use.	
Phone: 919-831-2240	Fax: 919-831-2240	email: info@immunoreagents.com

Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.

Trademark:

DyLight® is a trademark of Thermo Fisher Scientific, Inc. and its subsidiaries.