

## Product Specifications

Product:	LigaTrap™ Goat IgG Purification Resin
Description:	LigaTrap™ Goat IgG Purification Resin
Part Number:	LT-136
Concentration:	5 ml of LigaTrap™ Goat IgG Purification Resin
Amount:	5 ml
Form:	~50% Resin Slurry
Purification:	Binding capacity is greater than 15mgs of Goat IgG per ml resin
Buffer:	Not Applicable
Storage:	2-8 °C
Applications:	Affinity Purification
Disclaimer:	For <i>in vitro</i> Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. <b>Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.</b>