

## Product Specifications

|                |  |
|----------------|--|
| Product:       | LigaTrap™ Rat IgG Purification Resin   |
| Description:   | LigaTrap™ Rat IgG Purification Resin   |
| Part Number:   | LT-138   |
| Concentration: | 5 ml of LigaTrap™ Rat IgG Purification Resin   |
| Amount:        | 5 ml   |
| Form:          | ~50% Resin Slurry  |
| Purification:  | Binding capacity is greater than 20mgs of Rat IgG per ml resin   |
| Buffer:        | Not Applicable   |
| Storage:       | 2-8 °C   |
| Applications:  | Affinity Purification  |
| Disclaimer:    | For <i>in vitro</i> Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. <b>Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.</b> |