

Product Specifications

Product:	LigaTrap™ Chicken IgY Purification Resin
Description:	LigaTrap™ Chicken IgY Purification Resin
Part Number:	LT-142
Concentration:	5 ml of LigaTrap™ Chicken IgY Purification Resin
Amount:	5 ml
Form:	~50% Resin Slurry
Purification:	Binding capacity is greater than 15mgs of chicken IgY per ml resin
Buffer:	Not Applicable
Storage:	2-8 °C
Applications:	Affinity Purification
Disclaimer:	For <i>in vitro</i> Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.