

Product Specifications

Product:	LigaTrap™ Llama IgG Purification Resin
Description:	LigaTrap™ Llama IgG Purification Resin
Part Number:	LT-144
Concentration:	5 ml of LigaTrap™ Llama IgG Purification Resin
Amount:	5 ml
Form:	~50% Resin Slurry
Purification:	Binding capacity is greater than 15mgs of Llama IgG per ml resin. Binds IgG1, IgG2, and IgG3, isotypes
Buffer:	Not Applicable
Storage:	2-8 °C
Applications:	Affinity Purification
Disclaimer:	For <i>in vitro</i> Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.