

Product Specifications

Product:	Universal UltraPolymer Cocktail, Goat anti-Rabbit/Mouse IgG (H&L) – HRP
Description:	Universal UltraPolymer Cocktail, Goat anti-Rabbit/Mouse IgG (H&L) conjugated to HRP, affinity purified
Part Number:	UNIHRP-050
Concentration:	n/a
Amount:	50 ml
Conjugate:	Horseradish Peroxidase
Form:	Liquid
Purification:	Affinity purified using solid phase Rabbit IgG and Mouse IgG
Purity:	Affinity purified antibody is > 95% based on SDS-PAGE
Host:	Goat
Immunogen:	Purified Rabbit IgG, whole molecule; Purified Mouse IgG, whole molecule
Buffer:	PBS, 1% BSA & 0.1% proclin 150
Storage:	2-8 °C
Specificity:	Based on IEP, the Goat anti-Rabbit IgG immunopolymer antibody reacts with heavy (γ) chains on rabbit IgG, light chains on all rabbit immunoglobulins; the Goat anti-Mouse IgG immunopolymer antibody reacts with heavy (γ) chains on mouse IgG, light chains on
Cross Reactivity:	Based on IEP, no cross-reactivity is observed from Goat-anti-Rabbit IgG to non-immunoglobulin rabbit serum proteins, serum proteins from human, mouse or rat; no cross-reactivity is observed from Goat-anti-Mouse IgG to non-immunoglobulin mouse serum proteins, serum proteins from bovine, horse, human, pig or rabbit
Country of Origin:	Goat serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.
Applications:	For use as a standard or control in most immunoassay formats.
Disclaimer:	For <i>in vitro</i> Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.