

Product Data Sheet

Product:	PepMix Human (ACTL8)
Product code:	PM-ACTL8
Source:	Human Actin-like protein 8
Protein ID:	Q9H568
Gene ID:	ACTL8
No. of peptides:	89 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	366 aa
Sequence:	MAARTVIIDH GSGFLKAGTA GWNEPQMVPF NIVNYLPCKE NPGPSYARRR VSLGIDICHP DTFSYPIERG RILNWEQVQY LWSFVLENHR REQEVPPVII TETPLREPAD RKKMLEILFE LLHVPSVLLA DQLQMSLYAS GLLTGVVVDS GYGLTRVQPF HQGRPLPASG KTLFAGQDL SAYLLKSLFK EDCDRRCLFQ LETVAVTQMN KCYVPQNLGE ALDFRERQQS ALDESNITYQL PDGSRVELTP MQRVAPEMFF SPQVFEQPGP SIPRAIVESV ESCEISLRPL LVSHVMACGG NTLYPGFTKR LFRELMGDHV SSTKATVWEG SNRNFSVWLG ASVVAHLSTY QSEWMSREEY GEHMRM
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/
Use recommended within 6 months	