

Product Data Sheet

Product: PepMix Human (AFP)
Product code: PM-AFP
Source: Human
 Alpha-fetoprotein
Protein ID: P02771
Gene ID: AFP
No. of peptides: 150 peptides
Amount: 25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10⁸ cells)
Length: 609 aa

Sequence:

```

MKWVESIFLI FLLNFTEVRT LHRNEYGIAS ILDSYQCTAE ISLADLATIF FAQFVQEATY
KEVSKMVKDA LTAIEKPTGD EQSSGCLENQ LPAFLEELCH EKEILEKYGH SDCCSQSEEG
RHNCFLAHKK PTPASIPLFQ VPEPVTSCFA YEEDRETFMN KFIYEIARRH PFLYAPTILL
WAARYDKIIP SCCKAENAVE CFQTKAATVT KELRESSLLN QHACAVMKNF GTRTFQAITV
TKLSQKFTKV NFTEIQKLVV DVAHVHEHCC RGDVLDCLQD GEKIMSYICS QQDTLSNKIT
ECCKLTTLER GQCIIHAEND EKPEGLSPNL NRFLGDRDFN QFSSGEKNIF LASFVHEYSR
RHPQLAVSVI LRVAKGYQEL LEKCFQTENP LECQDKGEEE LQKYIQESQA LAKRSCGLFQ
KLGEYYLQNA FLVAYTKKAP QLTSSSELMAI TRKMAATAAT CCQLSEDKLL ACGEGAADII
IGHLCIRHEM TPVNPVGVCQ CTSSYANRRP CFSSLVVDTE YVPPAFSDDK FIFHKDLCQA
QGVALQTMKQ EFLINLVKQK PQITEEQLEA VIADFSGLLE KCCQGQEQEV CFAEEGQKLI
SKTRALGV
  
```

Purity: Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance: white to off-white lyophilisate
Remark: The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions: Store at - 20°C
Application protocol: <https://jpt.com/support-contact/resources/application-protocols/>

Use recommended within 6 months