

Product Data Sheet

Product:	PepMix BKV (small T antigen)
Product code:	PM-BKV-STA
Source:	BK Polyomavirus (BKV) Small T antigen
Protein ID:	P03082
No. of peptides:	41 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	172 aa
Sequence:	MDKVLNREES MELMDLLGLE RAAWGNLPLM RKAYLRKCKE FHPDKGGDED KMKRMNTLYK KMEQDVKVAH QPDFGTWSSS EVCADFPLCP DTLYCKEWPI CSKKPSVHCP CMLCQLRLRH LNRKFLRKEP LVWIDCYCID CFTQWFGLDL TEETLQWWVQ IIGETPFRDL KL
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months