

Product Data Sheet

Product:	PepMix BKV (capsid protein VP2)
Product code:	PM-BKV-VP2
Source:	BK Polyomavirus (BKV) Minor capsid protein - VP2
Protein ID:	P03094
No. of peptides:	85 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	351 aa
Sequence:	MGAALALLGD LVASVSEAAA ATGFSVAEIA AGEAAAAIEV QIASLATVEG ITSTSEAIAA IGLTPQTYAV IAGAPGAIAG FAALIQTVSG ISSLAQVGYR FFSDWDHKVS TVGLYQQSGM ALELFPNDEY YDILFPGVNT FVNNIQYLDP RHWGPSLFAT ISQALWHVIR DDIPSITSQE LQRRTEREFR DSLARFLEET TWTIVNAPIN FYNYIQQYYS DLSPIRPSMV RQVAEREGTR VHFGHTYSID DADSIEEVTQ RMDLRNQQSV HSGEFIEKTI APGGANQRTA PQWMLPLLLG LYGTVTPALE AYEDGPNQKK RRVSRGSSQK AKGTRASAKT TNKRRSRSSR S
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months