

Product Data Sheet

Product:	PepMix BKV (capsid protein VP2, isoform VP3)
Product code:	PM-BKV-VP3
Source:	BK Polyomavirus (BKV) Minor capsid protein - VP3
Protein ID:	P03094-2
No. of peptides:	56 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	232 aa
Sequence:	MALELFNPDE YYDILFPGVN Tfvnniqyld PRHWGPSLFA TISQALWHVI RDDIPSITSQ ELQRRTERFF RDSLARFLEE TTWTIVNAPI NFYNYIQYY SDLSPIRPSM VRQVAEREGT RVHFGHTYSI DDADSIEEVT QRMDLRNQQS VHSGEFIEKT IAPGGANQRT APQWMLPLLL GLYGTVTPAL EAYEDGPNQK KRRVSRGSSQ KAKGTRASAK TTNKRRSRSS RS
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months