

Product Data Sheet

Product:	PepMix Borrelia (OspA protein) Ultra
Label name:	Borrelia (OspA protein) Ultra
Product code:	PM-BORR-OSPA-1
Source:	Borrelia Outer surface protein A
Protein ID:	-
Gene ID:	ospA
No. of peptides:	200 peptides (Selected immunogenic peptides (15mers))
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	diverse
Sequence:	see additional product documentation
Purity:	Research Plus Grade: each peptide is guaranteed to be main product (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C.
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months