

Product Data Sheet

| | |
|------------------------------|--|
| Product: | PepMix Human (Brachyury) |
| Product code: | PM-BRAC |
| Source: | Human Transcription factors/activators |
| Protein ID: | O15178 |
| Gene ID: | T |
| No. of peptides: | 106 peptides |
| Amount: | 25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells) |
| Length: | 435 aa |
| Sequence: | MSSSPGTESAG KSLQYRVDHL LSAVENELQA GSEKGDPTER ELRVGLEESE LWLRFKELTN EMIVTKNGRR MFPVLKVNVS GLDPNAMYSF LLDFVAADNH RWKYVNGEWV PGGKPEPQAP SCVYIHPDSP NFGAHWMKAP VSFSKVKLTN KLNGGGQIML NSLHKYEPRI HIVRVGGPQR MITSHCFPET QFIAVTAYQN EEITALKIKY NPPAKAFLDA KERSDHKEMM EEPGDSQQPG YSQWGWLLPG TSTLCPANP HPQFGGALS PSTHSCDRYP TLRSHRSPY PSPYAHRNNS PTYSDNSPAC LSMLQSHDNW SSLGMPAHP MLPVSHNASP PTSSSQYPSL WSVSNGAVTP GSQAAAVSNG LGAQFFRGSP AHYTPLTHPV SAPSSSGSPL YEGAAAATDI VDSQYDAAAQ GRLIASWTFV SPFSM |
| Purity: | Development Grade: each peptide purified to > 70% (HPLC/MS) |
| Appearance: | white to off-white lyophilisate |
| Remark: | The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only. |
| Storage conditions: | Store at - 20°C |
| Application protocol: | https://jpt.com/support-contact/resources/application-protocols/ |

Use recommended within 6 months