

Product Data Sheet

Product:	PepMix Collection Breast Cancer
Product code:	PM-C-Brea
Source:	Human Selected proteins
Protein ID:	see datasheet
No. of peptides:	1018 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	diverse
Sequence:	8 different TAA PepMixes + 2 Controls (15mers, overlap 11): Human (CEA), Human (Claudin-6), Human (ErbB2_ECD), Human (ErbB2_ICD), Human (MAGEA3), Human (Mammaglobin A), Human (Mucin-1), Human (Stromelysin-3, MMP11), CEF Pool (extended), Human (Actin)
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months