

Product Data Sheet

Product:	PepMix Collection HCMV (2x1µg/peptide)
Product code:	PM-C-HCMV-1
Source:	Human Cytomegalovirus (HCMV) Selected proteins
Protein ID:	see datasheet
No. of peptides:	2743 peptides
Amount:	2 x aliquots of each peptide pool (1µg /peptide /pool)
Length:	diverse
Sequence:	19 different PepMixes, 2 aliquots of 19 peptide pools each containing 1µg per peptide /tube
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months