

Product Data Sheet

Product:	PepMix Candida (MP65)
Product code:	PM-CAND-MP65
Source:	Candida Albicans (Yeast) Mannoprotein MP65
Protein ID:	Q9HEP1
Gene ID:	mp65
No. of peptides:	92 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	379 aa
Sequence:	MLFKSFVTFV LANALAAPL AHQHHQHKEE KRAVHVVTIT NVVVVTIGNG DQTTTFAAPS VAADSSVSVS VNTEPPQNHP TTTQDVASAS TYPSSTDGSA ASSSAAASSS SQAGSEPSGG VGSGGAKGIT YSPYSDNGGC KSSSQIASEI AQLSGFNVIR LYGVDCDQVA AVLIAKTSSQ KIFAGIFDVS SITSGIESLA EAVKSICGSW DDIYTVSIGN ELVNAGSATP SQIKAYVEEG RKALKAAGYT GPVVSVDTFI AVINNPDLCD YSDYMAVNAH AFFDGHVAAE NSGAWVLQOI QRVWTACGGK KNVLITETGW PSRGDSNGVA VPSKSNQQA ISSIKSSCGA SAILFTAFND LWKADGPYNA EKYWGIYSN
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months