

Product Data Sheet

Product:	CEF Pool (extended)
Label name:	CEF Pool (extended)
Product code:	PM-CEF-E-1
Source:	Epstein-Barr Virus (EBV) Human Cytomegalovirus (HCMV) Influenza A Selected proteins
No. of peptides:	32 peptides (Defined HLA class I-restricted T-cell epitopes from Cytomegalovirus, Epstein-Barr virus, and Influenza virus.)
Amount:	15 nmol (approx. 25µg) / peptide for stimulation of $2,5 \times 10^8$ cells
Length:	diverse
Sequence:	see additional product documentation
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months