

## Product Data Sheet

<b>Product:</b>	CEF Pool (extended)
<b>Product code:</b>	PM-CEF-E-2
<b>Source:</b>	Epstein-Barr Virus (EBV) Human Cytomegalovirus (HCMV) Influenza A Selected proteins
<b>No. of peptides:</b>	32 peptides (Defined HLA class I-restricted T-cell epitopes from Cytomegalovirus, Epstein-Barr virus, and Influenza virus.)
<b>Amount:</b>	120 nmol /peptide
<b>Length:</b>	diverse
<b>Sequence:</b>	see additional product documentation
<b>Purity:</b>	Development Grade: each peptide purified to > 70% (HPLC/MS)
<b>Appearance:</b>	white to off-white lyophilisate
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

**Use recommended within 6 months**