

Product Data Sheet

| | |
|------------------------------|---|
| Product: | CEF Pool (extended) |
| Product code: | PM-CEF-E-3 |
| Source: | Epstein-Barr Virus (EBV) Human Cytomegalovirus (HCMV) Influenza A Selected proteins |
| No. of peptides: | 32 peptides (Defined HLA class I-restricted T-cell epitopes from Cytomegalovirus, Epstein-Barr virus, and Influenza virus.) |
| Amount: | 15 nmol (approx. 25µg) / peptide for stimulation of $2,5 \times 10^8$ cells |
| Length: | diverse |
| Sequence: | see additional product documentation |
| Purity: | Trial Grade: each peptide purified to > 90% (HPLC/MS) |
| Appearance: | white to off-white lyophilisate |
| Remark: | The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only. |
| Storage conditions: | Store at - 20°C |
| Application protocol: | https://jpt.com/support-contact/resources/application-protocols/ |

Use recommended within 6 months