

Product Data Sheet

Product: CEF Pool (extended)

Product code: PM-CEF-E-4

Source: Epstein-Barr Virus (EBV)

Human Cytomegalovirus (HCMV)

Influenza A
Selected proteins

No. of peptides: 32 peptides (Defined HLA class I-restricted T-cell epitopes from Cytomegalovirus, Epstein-Barr virus, and

Influenza virus.)

Amount: 120 nmol /peptide

Length: diverse

Sequence: see additional product documentation

Purity: Trial Grade: each peptide purified to > 90% (HPLC/MS)

Appearance: white to off-white lyophilisate

Remark: The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.

Storage conditions: Store at - 20°C

Application protocol: https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

