

Product Data Sheet

Product:	CEF Pool (standard)
Product code:	PM-CEF-S-4
Source:	Epstein-Barr Virus (EBV) Human Cytomegalovirus (HCMV) Influenza A Selected proteins
No. of peptides:	23 peptides (Defined HLA class I-restricted T-cell epitopes from Cytomegalovirus, Epstein-Barr virus, and Influenza virus.)
Amount:	120 nmol /peptide
Length:	diverse
Sequence:	see additional product documentation
Purity:	Trial Grade: each peptide purified to > 90% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months