

Product Data Sheet

Product:	CEFT Pool
Product code:	PM-CEFT-2
Source:	Clostridium Tetani Epstein-Barr Virus (EBV) Human Cytomegalovirus (HCMV) Influenza A Selected proteins
No. of peptides:	27 peptides (Defined HLA class I or II restricted T-cell epitopes from Cytomegalovirus, Epstein-Barr virus, Influenza virus or Clostridium tetani.)
Amount:	120 nmol /peptide
Length:	diverse
Sequence:	see additional product documentation
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months