

Product Data Sheet

Product:	CEFT MHC-II Pool
Product code:	PM-CEFT-MHC-II-3
Source:	Clostridium Tetani Epstein-Barr Virus (EBV) Human Cytomegalovirus (HCMV) Influenza A Selected proteins
No. of peptides:	14 peptides (Defined HLA class II restricted T-cell epitopes from Cytomegalovirus, Epstein-Barr virus, Influenza virus or Clostridium tetani.)
Amount:	15 nmol (approx. 25µg) / peptide for stimulation of 2,5 x 10 ⁸ cells
Length:	diverse
Sequence:	see additional product documentation
Purity:	Trial Grade: each peptide purified to > 90% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months