

Product Data Sheet

Source:

Product: CEFT MHC-II Pool
Product code: PM-CEFT-MHC-II-4

Clostridium Tetani Epstein-Barr Virus (EBV)

Human Cytomegalovirus (HCMV)

Influenza A Selected proteins

No. of peptides: 14 peptides (Defined HLA class II restricted T-cell epitopes from Cytomegalovirus, Epstein-Barr virus,

Influenza virus or Clostridium tetani.)

Amount: 120 nmol /peptide

Length: diverse

Sequence: see additional product documentation

Purity: Trial Grade: each peptide purified to > 90% (HPLC/MS)

Appearance: white to off-white lyophilisate

Remark: The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.

Storage conditions: Store at - 20°C

Application protocol: https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

