

Product Data Sheet

Product:	PepMix EBV (BMLF1)
Product code:	PM-EBV-BMLF1
Source:	Epstein-Barr Virus (EBV) strain B95-8 BMLF1 (mRNA export factor EB2)
Protein ID:	Q04360
Gene ID:	BMLF1
No. of peptides:	117 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	479 aa
Sequence:	MVPSQRLSRT SSISSNEDPA ESHILELEAV SDTNTDCDLD PMEGSEEHST DGEISSSEEE DEDPTPAHAI PARPSSVVIT PTSASFVIPR KKWDLQDKTV TLHRSPLCRD EDEKEETGNS SYTRGHKRRR GEVHGCTDES YGKRRHLPPG ARAPRAPRAP RVPRAPRSPR APRSNRATR PRSESRGAGR STRKQARQER SQRPLPNKPW FDMSLVKPVV KITFVTLPSL LASLTLEPIQ DPFLQSMMLAV AAHPEIGAWQ KVQPRHELRR SYKTLREFFT KSTNKDTWLD ARMQAIQNAG LCTLVAMLEE TIFWLQEITY HGDLPAPAE DILLACAMSL SKVILTKLKE LAPCFLPNTR DYNFVKQLFY ITCATARQNK VVETLSSSYV KQPLCLLAAY AAVAPAYINA NCRRRHDEVE FLGHYIKNYN PGTLSLLTE AVETHTRDCR SASCSRLVRA ILSPGTGSLG LFFVPGLNQ
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months