

Product Data Sheet

Product:	PepMix EBV (BMRF1)
Product code:	PM-EBV-BMRF1
Source:	Epstein-Barr Virus (EBV) strain B95-8 DNA polymerase processivity factor BMRF1
Protein ID:	P03191
Gene ID:	BMRF1
No. of peptides:	99 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	404 aa
Sequence:	METTQTLRFK TKALAVLSKC YDHAQTHLKG GVLQVNLLSV NYGGPRLAAV ANAGTAGLIS FEVSPDAVAE WQNHQSPEEA PAAVSFRNLA YGRTCVLGKE LFGSAVEQAS LQFYKRPQGG SRPEFVKLTM EYDDKVS KSH HTCALMPYMP PASDRLRNEQ MIGQVLLMPK TASSLQKWAR QQGSGGVKVT LNPDLVTTY TSGEACLTD YKPLSVGPYE AFTGPPVAKAQ DVGAVEAHVV CSVAADSLAA ALSLCRIPAV SVPILRFYRS GIIAVVAGLL TSAGDLPLDL SVILFNHASE EAAASTASEP EDKSPRVQPL GTGLQQRPRH TVSPSPSPPP PPRTPTWESP ARPETPSPAI PSHSSNTALE RPLAVQLARK RTSSEARQKQ KHPKKVKQAF NPLI
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months