

Product Data Sheet

Product: PepMix EBV (BZLF1)

Product code: PM-EBV-BZLF1

Source: Epstein-Barr Virus (EBV)

strain B95-8

BZLF1 Trans-activator protein

 Protein ID:
 P03206

 Gene ID:
 BZLF1

No. of peptides: 59 peptides (15mers with 11 aa overlap)

Amount: 25 μg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10⁸ cells)

Length: 245 aa

Sequence: MMDPNSTSED VKFTPDPYQV PFVQAFDQAT RVYQDLGGPS QAPLPCVLWP VLPEPLPQGQ

LTAYHVSTAP TGSWFSAPQP APENAYQAYA APQLFPVSDI TQNQQTNQAG GEAPQPGDNS TVQTAAAVVF ACPGANQGQQ LADIGVPQPA PVAAPARRTR KPQQPESLEE CDSELEIKRY KNRVASRKCR AKFKQLLQHY REVAAAKSSE NDRLRLLLKQ MCPSLDVDSI IPRTPDVLHE

DLLNF

Purity: Development Grade: each peptide purified to > 70% (HPLC/MS)

Appearance: white to off-white lyophilisate

Remark: The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.

Storage conditions: Store at - 20°C

Application protocol: https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

