

## Product Data Sheet

**Product:** PepMix EBV (EBNA3c)  
**Product code:** PM-EBV-EBNA3c  
**Source:** Epstein-Barr Virus (EBV)  
 strain AG876  
 Epstein-Barr nuclear antigen (EBNA)  
**Protein ID:** Q69140  
**Gene ID:** EBNA6  
**No. of peptides:** 265 peptides (15mers with 11 aa overlap)  
**Amount:** 25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10<sup>8</sup> cells)  
**Length:** 1069 aa

**Sequence:**

```

MESFEGEGDS IQSPDNARGD DVQNTGEHIQ DPGPGPSTGG ASEGLVQNEP DSRDQQSRGQ
RRGDENRGWM QRIRRRRRRR AALSGHLLDM EDNVPPWFPP HDITPYVARN IRDAACQAVK
HSHLQALSNL ILDSGLDTQH LLCFVMAARQ RLQDIRRGPL VVEGGVGWRH WLLTSPSRSW
SMGYRTATLR TLTPVPNRVG ADSIMLTATF GCQNGALAIN TFSATVWIPP PAGPREQERY
AREAEVRFLR GKWQRRFRRI FDLIELCGSL HHVWQNMLQT EENLLDFVRF MGVMSSCNSS
SVNYWFHKT I GNFKPYYPWN APPNENPYHA RRGIKEQVIQ KAFLKAQRQG LSMLATGGGP
RGDATSETSS DEDTGRQGS D VELESSDDEL PYIDPNMEPV QQRVPMFVSR VPVRKPRTL P
WPTPKTHPVK RTIVKTSYRS DEAEAAQSTP ERPGPSKQPS EPVEPAHTTP AGRSTVILHE
PPREPEAVSF KPPPPPSRRR RGACVVYDDD IIEVIDVETT EEETTSMQRQ PPLGQQPPPP
VISTGSAMSS SHTDPSVTQP SKPHRKPQDG FQRSGRRQKR AMPPPVSPSD AGPPSTRPRV
MAPPSTGPRV MATPSTGPRD MAPPSTGPRD MAPPSTGPRD MAPPSTGPRD MAPTVVHMFT
RERLLTQSTG PAPRSFWEMR AGRDAPKIQQ EPSSQQQPAT QSTPPCQSWV PSVYVLPADV
AGNAQPLQIS HLSSMSPTQP ISHEEQPRYE DPDTPLDLSL HPDTATLPPT QDLYPGREDL
QATQAPYPGY EEPRPPQAPF VGDYGFVQIP SAQWEPHPSQ GTYQGHIDPQ LPAALDLGPE
QPRFPQDPYV YSGGQLSSCP GYAGPWPSRP QHPRYRHTLA LWPREPRHGH SQGPWKPWSA
HLPPQWDGSA GHGQDQVSQF PHLHSETGPP RLQLSSVPPV LYPQPLVSSS APSWSSPQPR
APIRPIPTRF PPPPMPLQDS MAVGCDSSGT ACPSMPFASD YSQGAFTPLD INAPTPKSPR
VEESSHGPAR CSQATSEAQE ILSDNSEISV FPKDAKQTDY DASTESELD
  
```

**Purity:** Development Grade: each peptide purified to > 70% (HPLC/MS)  
**Appearance:** white to off-white lyophilisate  
**Remark:** The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.  
**Storage conditions:** Store at - 20°C  
**Application protocol:** <https://jpt.com/support-contact/resources/application-protocols/>

**Use recommended within 6 months**