

## **Product Data Sheet**

**Product:** EFX Ultra SuperStim Pool

Product code: PM-EFX-1

Source: Other/None

Selected proteins

No. of peptides: 146 peptides

**Amount:** 15 nmol (approx. 25μg) / peptide for stimulation of 2,5 x 10<sup>8</sup> cells

**Length:** diverse

**Sequence:** see additional product documentation

**Purity:** Development Grade: each peptide purified to > 70% (HPLC/MS)

**Appearance:** white to off-white lyophilisate

**Remark:** The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.

**Storage conditions:** Store at - 20°C

**Application protocol:** https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

