

## Product Data Sheet

**Product:** PepMix HAdV-3 (hexon protein)  
**Product code:** PM-HAdV3  
**Source:** Human Adenovirus (HAdV)  
 type B serotype 3  
 Hexon protein  
**Protein ID:** P36849  
**Gene ID:** L3  
**No. of peptides:** 234 peptides (15mers with 11 aa overlap)  
**Amount:** 25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10<sup>8</sup> cells)  
**Length:** 944 aa

**Sequence:**

```

MATPSMMPQW AYMHIAGQDA SGYLSPGLVQ FARATDTYFS MGNKFRNPTV APTHDVTTDR
SQRLMLRFVP VDREDNTYSY KVRVTLAVGD NRVLDMASTF FDIRGVLDGR PSFKPYSGTA
YNSLAPKGAP NTSQWIVTTN GDNAVTTTTN TFGIASMKGG NITKEGLQIG KDITTEGEE
KPIYADKTYQ PEPQVGEESW TDTDGTNEKF GGRALKPATN MKPCYGSFAR PTNIKGGQAK
NRKVKPTTEG GVETEEPDID MEFFDGRDAV AGALAPEIVL YTENVNLETP DSHVVYKPET
SNNSHANLGQ QAMPNRPNYI GFRDNFVGLM YNSTGNMGV LAGQASQLNA VVDLQDRNTE
LSYQLLLDSL GDRTRYFSMW NQAVDSYDPD VRIIENHGIE DELPNYCFPL NGIGPGHTYQ
GIKVKTDDTN GWEKDANVAP ANEITIGNNL AMEINIQANL WRSFLYSNVA LYLPDVYKYT
PPNITLPTNT NTYEYMNGRV VSPSLVDSYI NIGARWSLDP MDNVNPFNHH RNAGLRYSM
LLGNGRYVPF HIQVPQKFFA VKNLLLLPGS YTYEWNFRKD VNMVLQSSLG NDLRTDGATI
SFTSINLYAT FFPMAHNTAS TLEAMLRNDT NDQSFNDYLS AANMLYPIPA NATNIPISIP
SRNWAAFRGW SFTRLKTKET PSLGSGFDPY FVYSGSIPYL DGTFYLNHTF KKVAIMFDSS
VSWPGNDRLL SPNEFEIKRT VDGEYNVAQ CNMTKDWFLV QMLANYNIGY QGFYIPEGYK
DRMYSFFRNF QPMSRQVVDE VNYTDYKAVT LPYQHNSGF VGYLAPTMRQ GEPYPANYPY
PLIGTTAVKS VTQKKFLCDR TMWRIPFSSN FMSMGALTDL GQNMLYANSA HALDMTFEVD
PMDEPTLLYL LFEVFDVVRV HQPHRGVIEA VYLRTPFSSAG NATT
  
```

**Purity:** Development Grade: each peptide purified to > 70% (HPLC/MS)  
**Appearance:** white to off-white lyophilisate  
**Remark:** The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.  
**Storage conditions:** Store at - 20°C  
**Application protocol:** <https://jpt.com/support-contact/resources/application-protocols/>

**Use recommended within 6 months**