

Product Data Sheet

Product:	PepMix HBV (LEP) Ultra
Product code:	PM-HBV-LEPULTRA
Source:	Hepatitis B Virus (HBV) Large envelope protein
No. of peptides:	121 peptides (selected 15mers)
Amount:	25 μ g (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	diverse
Sequence:	see additional product documentation
Purity: Appearance: Remark:	Development Grade: each peptide purified to > 70% (HPLC/MS) white to off-white lyophilisate The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

JPT Peptide Technologies GmbH Volmerstrasse 5 12489 Berlin, Germany T +49-30-6392-5500 peptide@jpt.com | www.jpt.com

Geschäftsführer / Managing Directors Aurélien Claeyssen Sean Marett Dr. Sierk Pötting Sitz/Registered Office Berlin Gerichtsstand/Jurisdiction Berlin Amtsgericht Charlottenburg HRB 92692 B Ust.-IdNr./VAT-ID DE814044483 Bankverbindung/Bank Account Deutsche Bank Berlin Konto-Nr./Account-No. 915 165 500 IBAN DE84 1007 0000 0915 1655 00 BIC/SWIFT DEUTDEBBXXX

