

Product Data Sheet

Product: PepMix HEV (ORF3)

Product code: PM-HEV-ORF3

Source: Hepatitis E Virus (HEV)

Uncharacterized protein

Protein ID: E9N3B8

No. of peptides: 26 peptides (15mers with 11 aa overlap)

Amount: 25 μg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10⁸ cells)

Length: 113 aa

Sequence: MGSPCALGLF CCCSSCFCLC CPRHRPASRL AVVVGGAAAV PAVVSGVTGL ILSPSPSPIF

IQPTPSPPMS FHNPGLELAL GSRPAPLAPL GVXSPSAPPL PPAVDLPQLG LRR

Purity: Development Grade: each peptide purified to > 70% (HPLC/MS)

Appearance: white to off-white lyophilisate

Remark: The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.

Storage conditions: Store at - 20°C

Application protocol: https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

