

Product Data Sheet

| | |
|------------------------------|---|
| Product: | PepMix HIV-1 (Con B gag motif) |
| Product code: | PM-HIV-CONB |
| Source: | Human Immunodeficiency Virus (HIV) Con B gag motif |
| No. of peptides: | 123 peptides (15mers with 11 aa overlap) |
| Amount: | 25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells) |
| Length: | diverse |
| Sequence: | See Sequences under Product Documentation |
| Purity: | Development Grade: each peptide purified to > 70% (HPLC/MS) |
| Appearance: | white to off-white lyophilisate |
| Remark: | The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only. |
| Storage conditions: | Store at - 20°C |
| Application protocol: | https://jpt.com/support-contact/resources/application-protocols/ |

Use recommended within 6 months