

Product Data Sheet

Product:	PepMix Human (HMMR)
Product code:	PM-HMMR
Source:	Human Hyaluronan mediated motility receptor
Protein ID:	O75330
Gene ID:	HMMR
No. of peptides:	179 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	724 aa
Sequence:	MSFPKAPLKR FNDPSGCAPS PGAYDVKTLE VLKGPVSFQK SQRFKQQKES KQNLNVDKDT TLPASARKVK SSESKEQKN DKDLKILEKE IRVLLQERGA QDRRIQDLET ELEKMEARLN AALREKTSLS ANNATLEKQL IELTRTNELL KSKFSENGNQ KNLRIISLEL MKLRNKRETK MRGMMAKQEG MEMKLQVTQR SLEESQGKIA QLEGKLVSI EKEKIDEKSET EKLEYIEEI SCASDQVEKY KLDIAQLEEN LKEKNDEILS LKQSLEENIV ILSKQVEDLN VKCQLEKEK EDHVNRNREH NENLNAEMQN LKQKFILEQQ EREKLOQKEL QIDSLLOQEK ELSSSLHQKL CSFQEEMVKE KNLFEELKQ TLDELDKLQQ KEEQAERLVK QLEEEAKSRA EELKLEEKL KGKEAELEKS SAAHTQATLL LQEKYDSMVQ SLEDVTAQFE SYKALTASEI EDKLENSL QEKAAGKGN AEDVQHQILA TESSNQEYVR MLLDLQTKSA LKETEIKEIT VSFLQKITDL QNQLKQQEED FRKQLEDEEG RKAEKENTTA ELTEEINKWR LLYEELYNKT KPFQLQLDAF EVEKQALLNE HGAAQEQLNK IRDSYAKLLG HQNLKQKIKH VVKLKDENSQ LKSEVSKLRC QLAKKKQSET KLQEELNKVL GIKHFDPSKA FHESKENFA LKTPLKEGNT NCYRAPMECQ ESWK
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months