

Product Data Sheet

Product:	PepMix HPV16 (E5)
Product code:	PM-HPV16-E5
Source:	Human Papillomavirus (HPV) type 16 Probable protein E5
Protein ID:	P06927
Gene ID:	E5
No. of peptides:	18 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	83 aa
Sequence:	MTNLDTASTT LLACFLLCFC VLLCVCLLIR PLLLSVSTYT SLIILVLLW ITAASAFRCF IVYIIFVYIP LFLIHTHARF LIT
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months