

Product Data Sheet

Product: PepMix HPV 16 (E6)

Product code: PM-HPV16-E6

Source: Human Papillomavirus (HPV)

type 16

Protein E6

Protein ID: P03126

Gene ID: E6

No. of peptides: 37 peptides (15mers with 11 aa overlap)

Amount: 25 μg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10⁸ cells)

Length: 158 aa

Sequence: MHQKRTAMFQ DPQERPRKLP QLCTELQTTI HDIILECVYC KQQLLRREVY DFAFRDLCIV

YRDGNPYAVC DKCLKFYSKI SEYRHYCYSL YGTTLEQQYN KPLCDLLIRC INCQKPLCPE

EKQRHLDKKQ RFHNIRGRWT GRCMSCCRSS RTRRETQL

Purity: Development Grade: each peptide purified to > 70% (HPLC/MS)

Appearance: white to off-white lyophilisate

Remark: The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.

Storage conditions: Store at - 20°C

Application protocol: https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

