

Product Data Sheet

Product:	PepMix HPV 16 (L2)
Product code:	PM-HPV16-L2
Source:	Human Papillomavirus (HPV) type 18 Protein L2
Protein ID:	P03107
Gene ID:	L2
No. of peptides:	116 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	473 aa
Sequence:	MRHKRS AKRT KRASATQLYK TCKQAGTCPP DIIPKVEGKT IAEQILQYGS MGVFFGGLGI GTGSGTGGRT GYIPLGTRPP TATDTLAPVR PPLTVDPVGP SDPSIVSLVE ETSFIDAGAP TSVPSIPPDV SGFSITTSTD TTPAILDINN TVTTVTTHNN PTFTDPSVLQ PPTPAETGGH FTLSSSTIST HNYEEIPMDT FIVSTNPNTV TSSTPIPGSR PVARLGLYSR TTQQVKVVDV AFVTTPTKLI TYDNPAYEGI DVDNTLYFSS NDNSINIAPD PDFLDIVALH RPALTSRRTG IRYSRIGNKQ TLRTRSGKSI GAKVHYYDL STIDPAEEIE LQTITPSTYT TTSHAASPTS INNGLYDIYA DDFITDTSTT PVPSVPSTSL SGYIPANTTI PFGGAYNIPL VSGPDIPINI TDQAPSLIPI VPGSPQYTII ADAGDFYLHP SYYMLRKRK RLPYFFSDVS LAA
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months