

Product Data Sheet

Product: PepMix HPV18 (E5)

Product code: PM-HPV18-E5

Source: Human Papillomavirus (HPV)

type 18

Probable protein E5

Protein ID: P06792

Gene ID: E5

No. of peptides: 16 peptides (15mers with 11 aa overlap)

Amount: 25 μg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10⁸ cells)

Length: 73 aa

Sequence: MLSLIFLFCF CVCMYVCCHV PLLPSVCMCA YAWVLVFVYI VVITSPATAF TVYVFCFLLP

MLLLHIHAIL SLQ

Purity: Development Grade: each peptide purified to > 70% (HPLC/MS)

Appearance: white to off-white lyophilisate

Remark: The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.

Storage conditions: Store at - 20°C

Application protocol: https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

