

Product Data Sheet

Product:	PepMix HPV35 (E6)
Product code:	PM-HPV35-E6
Source:	Human Papillomavirus (HPV) type 35 Protein E6
Protein ID:	P27228
Gene ID:	E6
No. of peptides:	35 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	149 aa
Sequence:	MFQDPAERP Y KLHDL CNEVE ESIHEICLNC VYCKQELQRS EVYDFACYDL CIVYREGQPY GVCMKCLKFY SKISEYRWYR YSVYGETLEK QCNKQLCHLL IRCITCQKPL CPVEKQRHLE EKKR FHNIGG RWTGR CMSCW KPTRRETEV
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months