

Product Data Sheet

Product: PepMix HSV2 (E3 ICPO)
Product code: PM-HSV2-E3
Source: Herpes Simplex Virus
 strain HG52
 E3 ubiquitin-protein ligase ICPO
Protein ID: P28284
Gene ID: RL2
No. of peptides: 204 peptides (15mers with 11 aa overlap)
Amount: 25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10⁸ cells)
Length: 825 aa

Sequence:

```

MEPRPGTSSR ADPGPERPPR QTPGTQPAAP HAWGMLNDMQ WLASSDSEEE TEVGISDDDL
HRDSTSEAGS TDTEMFEAGL MDAATPPARP PAERQGSPTP ADAQGSCGGG PVGEEEEAEAG
GGGDVCAVCT DEIAPPLRCQ SFPCLHPFCI PCMKTWIPLR NTCPLCNTPV AYLIVGVTAS
GSFSTIPIVN DPRTRVEAEA AVRAGTAVDF IWTGNPRTAP RLSLSLGGHTV RALSPTPPWP
GTDEEDDDLA DVDYVPPAPR RAPRRGGGGA GATRGTSQPA ATRPAPPGAP RSSSSGGAPL
RAGVGS SGG GPAVA AVVPR VASLPPAAGG GRAQARRVGE DAAAAEGRTP PARQPRAAQE
PPIVISDSPP PSPRRPAGPG PLSFVSSSSA QVSSGPGGGG LPQSSGRAAR PRAAVAPRVR
SPPRAAAAPV VSASADAAGP APPAVPVDH RPRSRMTQA QTDTQAQSLG RAGATDARGS
GGPGAEGGPG VPRGTNTPGA APHAAEGAAA RPRKRRGSDS GPAASSSASS SAAPRSPLAP
QGVGAKRAAP RRAPDSDSGD RGHGPLAPAS AGAAPPSPSS SSQA AVAAAS SSSASSSSAS
SSSASSSSAS SSSASSSSAS SSSASSSAGG AGGSVASASG AGERRETSLG PRAAAPRGPR
KCARKTRHAE GGPEPGARDP APGLTRYLPI AGVSSVVALA PYVNKTVTGD CLPVLDMETG
HIGAYVVLVD QTGNVADLLR AAAPAWSRRT LLPEHARNCV RPPDYPTPPA SEWNSLWMTF
VGNMLFDQGT LVGALDFHGL RSRHPWSREQ GAPAPAGDAP AGHGE
  
```

Purity: Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance: white to off-white lyophilisate
Remark: The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions: Store at - 20°C
Application protocol: <https://jpt.com/support-contact/resources/application-protocols/>

Use recommended within 6 months