

Product Data Sheet

Product:	PepMix HSV2 (gD)
Product code:	PM-HSV2-gD
Source:	Herpes Simplex Virus strain HG52 Envelope protein
Protein ID:	Q69467
Gene ID:	gD
No. of peptides:	96 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	393 aa
Sequence:	MGRLTSGVGT AALLVVAVGL RVVCAKYALA DPSLKMADPN RFRGKNLPVL DQLTDPGGVK RVYHIQPSLE DPFQPPSIPI TVYYAVLERA CRSVLLHAPS EAPQIVRGAS DEARKHTYNL TIAWYRMGDN CAIPITVMEY TECPYNKSLG VCPIRTQPRW SYYDSFSAVS EDNLGFLMHA PAFETAGTYL RLVKINDWTE ITQFILEHRA RASCKYALPL RIPPAACLTS KAYQQGVTVD SIGMLPRFIP ENQRTVALYS LKIAGWHGPK PPYTSTLLPP ELSDTTNATQ PELVPEDPED SALLEDPAGT VSSQIPPNWH IPSIQDVAPH HAPAAPSNPG LIIGALAGST LAVLVIGGIA FWVRRRAQMA PKRLRLPHIR DDDAPPSHQP LFY
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months