

Product Data Sheet

Product:	PepMix Influenza A (HA/Hong Kong/1/1968 (H3N2))
Label name:	PepMix Inf. A (HA/HK1968/H3N2)
Product code:	PM-INFA-HAHK
Source:	Influenza A strain A/Hong Kong/1/1968 H3N2 Hemagglutinin
Protein ID:	Q91MA7
Gene ID:	HA
No. of peptides:	139 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	566 aa
Sequence:	MKTIIALSYI FCLALGQDLP GNDNSTATLC LGHHAVPNGT LVKTIITDDQI EVTNATELVQ SSSTGKICNN PHRILDGIDC TLIDALLGDP HCDVFQNETW DLFVERSKAF SNCYPYDVPD YASLRSLVAS SGTLEFITEG FTWTGVTQNG GSNACKRGP SGFFSRLNWL TKSGSTYPVL NVTMPNNDNF DKLYIWGVHH PSTNQEQTSL YVQASGRVTV STRRSQQTII PNIGSRPWVR GLSSRISIIYW TIVKPGDVLV INSNGLIAP RGYFKMRTGK SSIMRSDAPI DTCISECITP NGSIPNDKPF QNVNKITYGA CPKYVKQNTL KLATGMNRVVP EKQTRGLFGA IAGFIENGWE GMIDGWYGFR HQNSEGTGQA ADLKSTQAAI DQINGKLN RV IEKTNEKFHQ IEKEFSEVEG RIQDLEKYVE DTKIDLWSYN AELLVALENQ HTIDLTDSEM NKLFEKTRRQ LRENAEDMGN GCFKIYHKCD NACIESIRNG TYDHDVYRDE ALNNRFQIKG VELKSGYKDW ILWISFAISC FLLCVVLLGF IMWACQRGNI RCNICI
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months